

Connected
Biology
Everywhere



LABPAD[®] EVOLUTION



Use with:



**SMART[®]
COVID-19
TESTS**

Two meters in one handheld device

Unique patented lensless
technology.

Worry-free traceability


Test data (batch number, expiry date
and calibration data) automatically
read out from datamatrix.

Patient Identification

Connect to your smartphone,
PC or optional barcode reader.

Immediate setup with LabPad+

Complementary patient identification
APP for ANDROID and IOS



SMART[®] INR



avalun.com

Evoluteve panel of tests

CURRENTLY AVAILABLE

- SARS-COV2 Antigen Rapid Test
- TP/INR test

UNDER DEVELOPMENT:

- COVID-19 semi-quantitative IgG/IgM (Q1 2021)
- Prealbumin for undernutrition
- Cardiac markers

Designed to fit your practice



PATIENT
AT HOME



PRIMARY
CARE




LABORATORY



NURSING
HOMES

Mobile-Health ready



Patient identification



Test performance



Result broadcast

Onsite Near Patient Testing

Cloud services



Result processing

Centralized Clinical Information System


Advanced technology geared towards simplicity ...

All results in just three steps

K SMART®



T SMART®



... without compromise



Versatility

Two technologies in one single device allowing for a unique evolutive panel of tests.



Miniaturization

Handheld device to assist healthcare professionals in direct contact with patients wherever they are.



Reliability

Built in automatic quality controls to ensure the integrity of each result.

Features

Consumables	Ksmart® and Tsmart® Touchless ejection preventing for contamination
Memory	1 000 results
Connectivity	Bluetooth Low Energy USB B port
Power input	Built-in battery with charger
Operating conditions	15-32 °C 59-90 °F
Size	L 17.2 cm W 7.4 cm H from 2.0 to 4.6 cm
Weight	270 g



Lapad® Evolution box

Accessories (optional)

Bluetooth barcode reader
Bluetooth printer



Ordering information

CATALOG REFERENCE NUMBER:

LabPad® Evolution: 03760097080015
(with battery charger, carrying bag and user manual)
Ksmart® SARS-COV2 Antigen Rapid Test: 03760097080046
Tsmart® INR: box of 12 microcuvettes 03760097080121
box of 24 microcuvettes 03760097080138
box of 48 microcuvettes 03760097080145



AVALUN SAS

7, Parvis Louis Néel, 38000 Grenoble
France

contact@avalun.com



Find all the information on

avalun.com