

Connected  
Biology  
Everywhere



# LABPAD<sup>®</sup> EVOLUTION



Use with:



**SMART<sup>®</sup>  
COVID-19  
TESTS**

## Two meters in one handheld device

Unique patented lensless  
technology.

## Worry-free traceability


Test data (batch number, expiry date  
and calibration data) automatically  
read out from datamatrix.

## Patient Identification

Connect to your smartphone,  
PC or optional barcode reader.

## Immediate setup with LabPad+

Complementary patient identification  
APP for ANDROID and IOS



**SMART<sup>®</sup> INR**



[avalun.com](http://avalun.com)

# Evolutionary panel of tests

## CURRENTLY AVAILABLE

- SARS-COV2 Antigen Rapid Test
- TP/INR test

## UNDER DEVELOPMENT:

- COVID-19 semi-quantitative IgG/IgM (Q1 2021)
- Prealbumin for undernutrition
- Cardiac markers

## Designed to fit your practice



PATIENT AT HOME



PRIMARY CARE



LABORATORY



NURSING HOMES

## Mobile-Health ready



Patient identification

Test performance

Result broadcast

Onsite Near Patient Testing

## Cloud services



And/Or




Result processing

Centralized Clinical Information System


# Advanced technology geared towards simplicity ...

All results in just three steps

**K SMART®**



**T SMART®**



## ... without compromise



### Versatility

Two technologies in one single device allowing for a unique evolutive panel of tests.



### Miniaturization

Handheld device to assist healthcare professionals in direct contact with patients wherever they are.



### Reliability

Built in automatic quality controls to ensure the integrity of each result.

## Features

<b>Consumables</b>	Ksmart® and Tsmart® Touchless ejection preventing for contamination
<b>Memory</b>	1 000 results
<b>Connectivity</b>	Bluetooth Low Energy USB B port
<b>Power input</b>	Built-in battery with charger
<b>Operating conditions</b>	15-32 °C   59-90 °F
<b>Size</b>	L 17.2 cm W 7.4 cm H from 2.0 to 4.6 cm
<b>Weight</b>	270 g



Lapad® Evolution box

## Accessories (optional)

Bluetooth barcode reader  
Bluetooth printer



## Ordering information

### CATALOG REFERENCE NUMBER:

LabPad® Evolution: 03760097080015  
(with battery charger, carrying bag and user manual)  
Ksmart® SARS-COV2 Antigen Rapid Test: 03760097080046  
Tsmart® INR: box of 12 microcuvettes 03760097080121  
box of 24 microcuvettes 03760097080138  
box of 48 microcuvettes 03760097080145



**AVALUN SAS**

7, Parvis Louis Néel, 38000 Grenoble  
France

[contact@avalun.com](mailto:contact@avalun.com)



Find all the information on  
[avalun.com](http://avalun.com)